- 1. Unzip the "Information for GX-Force Leak Check Mode Upgrade" folder.
- 2. Navigate to the "USB C Program Loader" folder in the unzipped folder and double click the setup.exe file. Follow installation prompts.
- 3. Navigate to the "Setup Program" folder in the unzipped folder and double click the setup.exe file. Follow installation prompts.
- 4. Turn the GX-Force off if it is on.
- 5. Connect the GX-Force to the computer using the USB cable.



- 6. Launch the GX-Force Program Loader program from the computer's Start menu.
- 7. Turn on the GX-Force.

8. Click the Open button in the Program Loader interface.

RIKEN GX-Force Program Loader (Pno.06805)	1. <del></del> - 1.		×
PORT COM11 ~ Open Close	Main	🔿 Sub	
FILE			
SELECT			
STATUS			
			_
		GC	)
Ready			

9. Click the "Main" radio dot and then click "Select".

RIKEN GX-Force Program Loader (Pno.06805)	1. <del></del>		<
PORT COM11 V Open Glose	Main	🔿 Sub	
FILE			
SELECT			
STATUS			
			_
			-
		GO	
Ready			

10. Navigate to the Instrument Firmware\Main folder in the unzipped folder and double click the .mot file there.

9 Open				×
🔶 🐳 🕆 📙 « GX-Force > Instrument Firmware > Main	~	Ō	Search Main	Q
Organize 👻 New folder				
Videos ^ Name	Date modifie	d	Туре	Size
<ul> <li>OS (C:)</li> <li>D7554.mot</li> <li>HD2 (D:)</li> <li>USB20FD (E:)</li> <li>ISO 9000 (I:)</li> <li>Konica Scans (K:</li> <li>Engineering (N:)</li> <li>Departments (P:</li> <li>Users (V:)</li> <li>USB20FD (E:)</li> </ul>	2/7/2024 8:4:	AM	MOT File	98;
V <				>
File name:		~	S format(*.mot)	~
			Open	Cancel

11. Click Go. The program screen indicates the transfer is starting and the GX-Force screen says "UPDATE".

RIKEN G	X-Force Program	m <mark>Loader (</mark> F	no.06805)		. <del></del>		×
PORT	COM11 🖂	Open	Close	۲	Main	O Su	ub
FILE							
SELECT	N:\Eng Docu	ments\Unit s	oftwares\GX-Force\lr	nstrument Firr	mware\/	Main \075	54.mot
RC,R,							-
						C	GO
ON Line							

12. The program will indicate the transfer is complete. Click OK. The GX-Force will be turned off.



13. Click the "Sub" radio dot and then click "Select".

RIKEN GX-Force Program Loader (Pno.06805)	- \	□ ×
PORT COM11 V Open Close	🔘 Main	Sub
FILE		
SELECT		
STATUS		
		GO
Ready		

14. Navigate to the Instrument Firmware\Sub folder in the unzipped folder and double click the .mot file there.

9 Open			×
← → × ↑ 🦲 « GX-Force > Instrument Firmware > Sub	ٽ ~	Search Sub	م
Organize 👻 New folder			
Music Name	Date modified	Туре	Size
Videos	2/7/2024 8:49 AM	MOT File	40
<ul> <li>OS (C:)</li> <li>HD2 (D:)</li> </ul>			
■ USB20FD (E:) ■ ISO 9000 (I:)			
Engineering (N:)			
➡ Departments (P:) ➡ Users (V:)			
■ USB20FD (E:) 🗸 <			>
File name:	~	S format(*.mot) Open	∨ Cancel

15. Turn the GX-Force back on.

## 16. Click the Open button.

RIKEN GX-Force Program Loader (Pno.06805)	10		Х
PORT COM11 V Open Close	O Main	🖲 Su	ıb
FILE			
SELECT N:\Eng Documents\Unit softwares\GX-Force\Instr.	ument Firmware\S	Sub\0755	5.mot
STATUS			
			~

17. Click "Go". The program screen indicates the transfer is starting and the GX-Force screen says "UPDATE".

RIKEN G	GX-Force Progra	m Loader (P	no.06805)	1	
PORT	COM11 ~	Open	Close	🔿 Main	Sub
FILE					
SELECT	N:\Eng Docu	iments\Unit so	oftwares\GX-Force\l	nstrument Firmware\S	Gub\07555.m
STATUS					
RC,R,					
					GO
					40

18. The program will indicate the transfer is complete. Click OK. The GX-Force will be turned off.

